

Top End® Excelerator™ Stock Handcycle



| PROVIDER / CUSTOMER INFORMATION: | | |
|--|------------------|--|
| Provider Information: | | End User Information: |
| Date: | Contact Name: | Athlete Name: |
| Acct #: | Ship to Address: | |
| Contact: | City: | State: |
| Phone #: | Zip/Postal Code: | Phone: |
| Fax/Email: | Email: | |
| <input type="radio"/> Ship to Provider | | <input type="radio"/> Ship to End User |

| ORDER REQUEST TYPE: | | |
|-----------------------------|-----------------------------|---------------------------|
| <input type="radio"/> Quote | <input type="radio"/> Order | PO # <input type="text"/> |

| BASE MODEL: | | |
|--|-----------------------------|---------|
| <input checked="" type="radio"/> 10006 | Excelerator Stock Handcycle | \$4,499 |

| Assembly: <i>(see definitions)</i> | | |
|--|-----------|-----|
| <input checked="" type="radio"/> 10048 | Assembled | \$0 |

| SEAT WIDTH: | | |
|--|-----------------------------------|-----|
| <input checked="" type="radio"/> 15021 | 18" Seat Width - Royal Blue Frame | \$0 |

| CAMBER ANGLE: | | |
|--|-----------------|-----|
| <input checked="" type="radio"/> 45002 | 3 Degree Camber | \$0 |

| WHEELBASE: <i>(Front to Rear Axle)</i> | | |
|--|---------------|-----|
| <input checked="" type="radio"/> 75085 | 40" Wheelbase | \$0 |

| WHEEL SETS: | | |
|--|---|-----|
| <input checked="" type="radio"/> 50121 | 24" Front Wheels & 24" Rear Wheels with Quick Release Axles | \$0 |

| TIRES: | | |
|--|---------------------------|-----|
| <input checked="" type="radio"/> 50213 | 24" X 1.75" Treaded Tires | \$0 |

| FOOTREST: | | |
|--|-------------------------------|-----|
| <input checked="" type="radio"/> 60002 | Dual, Height/Angle Adjustable | \$0 |

| DRIVE CONTROL: | | |
|--|---|-----|
| CRANKS: | | |
| <input checked="" type="radio"/> 65054 | Crank Arms: 165mm Long by 12" Wide | \$0 |
| HAND PEDALS: | | |
| <input checked="" type="radio"/> 65137 | Medium 3D Printed with Power Plates | \$0 |
| SHIFTING: | | |
| <input checked="" type="radio"/> 65060 | Twist Activated Shifter with Internal Gearing | \$0 |
| BRAKING: | | |
| <input checked="" type="radio"/> 65080 | Reversing Activated Drum Brake | \$0 |
| <input checked="" type="radio"/> 65090 | Locking Parking Brake | \$0 |
| CHAINRING: | | |
| <input checked="" type="radio"/> 65100 | Single 38 Tooth Chainring | \$0 |

Top End® Excelerator™ Stock Handcycle



| UPHOLSTERY: | | |
|----------------------------------|--|-----|
| MATERIAL & COLOR: | | |
| <input checked="" type="radio"/> | 30001 Nylon, Black Upholstery | \$0 |
| SEAT UPHOLSTERY: | | |
| <input checked="" type="radio"/> | 30014 Adjustable Tension Seat Upholstery | \$0 |
| <input checked="" type="radio"/> | 25060 Hook and Loop Positioning Strap | \$0 |
| BACK UPHOLSTERY: | | |
| <input checked="" type="radio"/> | 30029 Fold - Over Back Upholstery | \$0 |
| SEAT CUSHION: | | |
| <input checked="" type="radio"/> | 25063 Seat Cushion | \$0 |

| UPGRADABLE OPTIONS AND ACCESSORIES: | | |
|--|-----------------------------|-----|
| <i>(Accessories will not be installed and may be shipped separately - See Definitions)</i> | | |
| HYDRATION: | | |
| <input checked="" type="radio"/> | 25012 Water Bottle and Cage | \$0 |

| STANDARD ACCESSORIES: | | |
|----------------------------------|-------------------|-----|
| <input checked="" type="radio"/> | 25061 Safety Flag | \$0 |
| <input checked="" type="radio"/> | 25062 Chain Guard | \$0 |

| HANDCYCLE SPECIFICATIONS: | |
|-------------------------------|--|
| Weight: | 50 pounds |
| Turning Radius: | 12 feet for 360 degree circle, 8 feet on a 3-point turn |
| Transfer Height: | 20" from ground with 2" cushion |
| Frame Style: | One piece main frame, upright position, 1-3/4" X .065" ovalized carbon steel tubing |
| Fork: | Carbon Steel |
| Propulsion / Shifting: | Internal front hub, twist activated shifter mounted on crank assembly |
| Cranks: | 165mm (length) X 12" (width), adjustable height |
| Back: | Adjustable height from 14"-18" |
| Footrest: | Dual adjustable height and angle, includes straps |
| Hand Pedals: | Ovalized ergonomic with frictionless bearings |
| Brake: | Reverse Activated Drum and Locking Parking Brake mounted on fork |
| Wheelbase: | 40" Axle to Axle |
| Camber: | 3 degree |
| Wheels: | 24" metal spoked wheels |
| Axles: | Quick release |
| Tires / Tubes: | 24" X 1.75" treaded tires (40 PSI) |
| Positioning Straps: | Hook and loop seat restraint included |
| Cushions: | Seat cushion included |
| Safety / Visibility: | Flag and wheel reflectors included |
| Chain Guard: | Full length, shock resistant composite material except when Mountain Drive is chosen. With Mountain Drive selection a half moon shaped black composite chainguard is used. |

| DEFINITIONS: | |
|------------------|--|
| Shipping: | Arrives factory assembled, rear wheels and some accessories will need to be installed. Shipping rate is \$325 for all shipments in the continental US. Freight quotes for shipments outside the continental US are available upon request. |

Warranty Info:

Top End warrants the frames when purchased new and unused to be free from defects in materials and workmanship for a period of three (3) years from the date of purchase from Top End or a dealer, with a copy of the seller's invoice required for coverage under this warranty. Top End warrants the upholstered materials (seat and back) and remaining components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of thirteen (13) months from date of purchase from Top End or a dealer, with a copy of the seller's invoice required for coverage under this warranty.